



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Spare paper holder, polished steel Q10260
<b>Product description</b>	VVS:776591240 NRF:6190933 RSK:8744046 EAN:5708590352831  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)