



PRESSALIT®

Product Item number	Wall-mounted toilet 27 4/8" R2099000
Color	VVS:612928000 NRF:6191026 RSK:7822154 EAN:5708590345284 HMI: 105700 White
Product description	Toilet, wall-mounted, distance from wall to the front: 27 4/8" Excl. mounting kit. For mounting directly on wall, a mounting kit from V&B must be ordered from Pressalit at an additional cost. For mounting on built-in cistern or Pressalit toilet lifter, mounting kit is included in the delivery. Toilet, wall-mounted, distance from wall to the front: 27 4/8". Excl. mounting kit.
Product capacity	Max load: 880 lbs.
Materials	Sanitary porcelain
Ergonomy	27 4/8" is recommended for wheelchair users as it provides sufficient space for a caregiver. Wall mounted toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a floor lift, and the caregiver's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Wall-mounted toilet 27 4/8"
Item number	R2099000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026