




Product Item number	Pressalit Choice Towel hook, curve, 2 pcs., polished steel Q10160
Product description	<p>VVS:775295340 NRF:6190929 RSK:8754883 EAN:5708590336879</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Use the towel hook in the bathroom and in the kitchen. That way, you can let the style be consistent throughout the house.</p>
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty	A 2-year warranty is provided for material and manufacturing defects.
Cleaning	<p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>