




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 950 ml, incl. stand, brushed steel PY30735
Product description	VVS: 774565815 NRF: 1438618 RSK: 8920392 EAN: 5708590392349  Disinfection dispenser in brushed stainless steel. Mounted on floor standing stand. Disinfection pump with sensor and 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Incl. mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com