



PRESSALIT®

Pressalit Style D2 1206

**Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Incl. design emblems, chrome and black**

Features

Lift-off, Soft close

Item number

1206111-DN2999
VVS:615099111
NRF:1476358
RSK:7824134
EAN:5708590397481

Pressalit Style D2 is designed for those who want more than just a toilet seat. Its slim, refined profile creates a light and exclusive expression that complements any contemporary bathroom.

With interchangeable side design emblems, you can create a cohesive look with the rest of your Pressalit Style bathroom accessories and tailor the details to match your personal taste. The seat comes with two sets of design emblems, making it easy to change the colour and appearance whenever you want a refresh. Additional design emblems in other colours and finishes are also available for purchase.

Pressalit Style D2 is available in six beautiful colours: white, matt white, black, matt black, matt dusty green and matt taupe. These shades allow you to create exactly the atmosphere you want in your bathroom.

Behind the design lies Pressalit's renowned quality, a 10-year warranty, and features such as Soft Close for a quiet, luxurious closing and Lift-off for easy and efficient cleaning. In addition, Pressalit SteadyFix ensures exceptional stability.

**Colour**

Black (RAL 9005)

Description




PRESSALIT toilet seat with cover, model "Pressalit Style D2", art. no. 1206 made of colour ingrained duroplast with antibacterial effect, with built-in soft close function, incl. DN2 universal hinge with lift-off in stainless steel.
- centre distance: 144-224 mm
- load - ring seat 240 kg
- 10 year guarantee



PRESSALIT®

Pressalit Style D2 1206

Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel
Incl. design emblems, chrome and black



Design type	Sandwich
Designer	Pressalit
Design year	2025
Materials	TOILET SEAT: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde). UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing ($< 0.01 \text{ mg/m}^3$) to emit significantly below the permitted EU limit of 0.08 mg/m^3. ANTIBACTERIAL: This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth. BUFFERS: EVA (copolymer made of ethylene and vinyl acetate). DAMPER KIT: Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic. MOUNT/HINGES: Stainless steel and thermoplastic.
Product capacity	Load - ring seat 240 kg
Weight	2,665 kg
Carton size	457x389x68 mm
Gross weight	3,089 kg
Crate size	469x401x348 mm



PRESSALIT®

Pressalit Style D2 1206

**Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

Incl. design emblems, chrome and black**

Weight per crate 15,445 kg

Qty per crate 5

Crates per pallet 18

Qty per pallet 90

Cleaning Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.

Warranty Your Pressalit toilet seat is subject to a 10-year product warranty against presentation of receipt. You should contact the dealer or store where you purchased the seat in the event of a complaint. In the event of complaints about defective materials and/or manufacturing defects and/or functional limitations ascertained within a 10-year period, following proper use and cleaning of the toilet seat (we reserve the right to examine the toilet seat at our factory), Pressalit will deliver spare parts or replace the toilet seat by an identical or similar model. If the color is no longer available, Pressalit reserves its right to replace the toilet seat by a similar model in the color white. Upon replacement or delivery of a spare part, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit may instead refund the purchase price.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com
