




PRESSALIT®

Product Item number	Pressalit Choice Towel hook, 2 pcs., matt white Q60100
Product description	<p>VVS:775295130 NRF:7001068 RSK:8754993 EAN:5708590370422 DB: 2173260</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom</p> <p>Use the towel hook in the bathroom and in the kitchen. That way, you can let the style be consistent throughout the house.</p>
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®