




<b>Product</b> <b>Item number</b>	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel PY30725
<b>Product description</b>	<p>VVS: 774565615          NRF: 1438614          RSK: 8920390          EAN: 5708590392745</p>  <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b> <b>WEEE</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>



# PRESSALIT®

<b>Product</b>	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel
<b>Item number</b>	PY30725

---

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®

8. juli 2026