



PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 large roll of toilet paper, polished steel PX10225
Product description	VVS: 776130140 NRF: 1438651 RSK: 8920360 EAN: 5708590392080 Toilet paper holder in large steel housing with space for rolls of max. Ø 280 x 100 mm. For wall mounting. The holder is lockable and includes mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com