




PRESSALIT®

Product Item number	Freestanding column 945 mm, for fixed height PLUS support arm R9430
Product description	VVS:615851990 NRF:6190667 EAN:5708590320427 NBS:N13:436  Free-standing column, grey, for fixed height PLUS support arm. Incl. mounting kit for floor. Free-standing column, grey, for fixed height PLUS support arm. Incl. mounting kit for floor.
Materials	Basic material: Steel Surface material: Rilsan
Features	Used with freestanding toilets and where it is not possible to wallmount a support arm. For fixed height support arm.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

pressalit.com

PRESSALIT®