



<b>Product</b>	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, matt black
<b>Item number</b>	PY40725
<b>Product description</b>	<p>VVS: 774565631 NRF: 1438615 RSK: 8920399 EAN: 5708590392769</p> <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b>	A 2-year warranty is provided for material and manufacturing defects.
<b>Cleaning</b>	<p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>