



PRESSALIT®

Product Item number	Pressalit Public Universal Disinfection dispenser with sensor, 950 ml, brushed steel PZ30732
Product description	VVS: 774566615 NRF: 1438660 RSK: 8920592 EAN: 5708590396361 Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning WEEE	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive. Subject to changes in the product range, technical changes and printing errors.

pressalit.com