



United Kingdom Conformity Assessment

United Kingdom RP:
 Pressalit Ltd.
 100, Longwater Avenue
 Green Park
 Reading RG2 6GP
 United Kingdom
 T: 01180923 2696
 uk@pressalit.com

This product is produced in compliance with all relevant requirements of:

- Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
- Council Directive 2006/42/EC on machinery
- Council Directive 2011/65/EC on RoHS

For the UK:

- Medical Device Regulation (MDR) 2002 on medical devices
- Supply of Machinery (Safety) Regulations 2008
- Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019
 ISO 13485:2016

The products are manufactured in accordance with the following standards:

IEC 60601-1:2005
 ISO 17966:2016

Ry, 28. juli 2026

Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S
 Responsible

Manufacturer: Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

SELECT TL3



R9123112

EAN (GTIN): 5708590372808

Nomenclature (GMDN): 63055

Nomenclature (EMDN): Y091299

Name

EN SELECT TL3 toilet bracket, for wall outlet, anthracite single molded cover. Electrically height adjustable with control panel (on support arm)

Intended purpose

EN Toilet lifters ensure optimum ergonomics at a toilet, for both seated and standing users with disabilities. Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.