



# PRESSALIT®

|                                |   |
|--------------------------------|---|
| <b>Product<br/>Item number</b> | Pressalit Public Classic Soap dispenser with sensor, 850 ml, polished steel<br>PY10715  |
| <b>Product description</b>     | VVS:774565240<br>NRF:1438641<br>RSK:8920382<br>EAN:5708590392264<br><br>Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.   |
| <b>Materials</b>               | European stainless steel (AISI 316L)  |
| <b>Warranty<br/>Cleaning</b>   | A 2-year warranty is provided for material and manufacturing defects.<br><br>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.<br><br>Subject to changes in the product range, technical changes and printing errors. |

[pressalit.com](http://pressalit.com)