



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Universal Soap dispenser, 1100 ml, brushed steel PZ30710
<b>Product description</b>	VVS:774566115 NRF:1438657 RSK:8920589 EAN:5708590396330  Stainless steel soap dispenser for wall mounting. One-handed, drip-free soap pump with 1100 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)