




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, polished steel PY10725
Product description	VVS: 774565640 NRF: 1438616 RSK: 8920383 EAN: 5708590392271  Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com