




PRESSALIT®

Fixed hinge

Item number

DA7999
VVS: 616810879
EAN: 5708590315041



Product capacity

Weight	0,148 kg
Gross weight	0,148 kg
Weight per crate	0 kg
Qty per crate	0
Qty per pallet	0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®