




PRESSALIT®

Product Item number	Pressalit Choice Towel holder, double, 60 x 12 cm, brushed brass Q70550
Product description	VVS:775711205 NRF:7001581 RSK:8920051 EAN:5708590381893 DB: 2425192  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com