




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Premium Toilet paper holder for 1 roll of toilet paper + 1 spare roll, matte black PX40210
<b>Product description</b>	VVS: 776130231 NRF: 1438647 RSK: 8920374 EAN: 5708590392691  <p>Toilet paper holder in steel housing with space for spare roll. For wall mounting. The holder is lockable and includes mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 304)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)