



PRESSALIT®

Product Item number	Pressalit Public Classic Waste bin 34 l, polished steel PY10910
Product description	VVS:777744340 NRF:1438610 RSK:8920385 EAN:5708590392295 Waste bin made of robust stainless steel. Stainless steel body with rounded edges. Equipped with removable, washable polypropylene fabric bag. Capacity approx. 34 liters. Installation can be done with glue or with two screws. There is also the option of mounting with extra anti-theft protection in the back wall. Incl. mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com