



PRESSALIT®

| | |
|--------------------------------------|---|
| Product Item number | Pressalit Public Premium Hygienic waste bin 4 l, matt black PX40925 |
| Product description | VVS:777743431 NRF:1438635 RSK:8920379 EAN:5708590392233 Wall-mounted sanitary bin with bag dispenser integrated on the side, removable 4 liter sanitary bin and hinged lid. Includes mounting kit and a box with plastic sanitary bags. |
| Materials | European stainless steel (AISI 304) |
| Warranty Cleaning | A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors. |

pressalit.com