




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser, 950 ml, brushed steel PY30730
Product description	<p>VVS:774565515 NRF:1438617 RSK:8920391 EAN:5708590392332</p>  <p>Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com