




PRESSALIT®

Product	Headrest for MSCT 1 and MSCT 4
Item number	R844521302
	VVS: 777542892 EAN: 5708590374932 HMI: 138482
Colour	Pewter grey
Product description	Headrest for MSCT 1 and MSCT 4 shower change trolley, 660 x 220 x 100 mm. Headrest, 660 x 220 x 100 mm, for MSCT 1 and MSCT 4 shower change trolley.
	

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®