




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, matt black PY40725
Product description	VVS: 774565631 NRF: 1438615 RSK: 8920399 EAN: 5708590392769  Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com