




<b>Product</b> <b>Item number</b>	Pressalit Public Classic Disinfection dispenser, 950 ml, polished steel PY10730
<b>Product description</b>	<p>VVS: 774565540 NRF: 1438620 RSK: 8920384 EAN: 5708590392288</p>  <p>Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b>	A 2-year warranty is provided for material and manufacturing defects.
<b>Cleaning</b>	<p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>