




# PRESSALIT®

|                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|---------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b>              | Pressalit Public Universal Disinfection dispenser with sensor, 950 ml, brushed steel<br>PZ30732                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product description</b>                        | <p>           VVS: 774566615<br/>           NRF: 1438660<br/>           RSK: 8920592<br/>           EAN: 5708590396361         </p>  <p>           Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.         </p> <p>           Inspection opening indicates filling level.         </p>                                     |
| <b>Materials</b>                                  | European stainless steel (AISI 316L)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Warranty</b><br><b>Cleaning</b><br><b>WEEE</b> | <p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p> |



# PRESSALIT®

|                    |                                                                                      |
|--------------------|--------------------------------------------------------------------------------------|
| <b>Product</b>     | Pressalit Public Universal Disinfection dispenser with sensor, 950 ml, brushed steel |
| <b>Item number</b> | PZ30732                                                                              |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®

8. juli 2026