



PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser, 950 ml, matt black PY40730
Product description	VVS:774565531 NRF:1438619 RSK:8920400 EAN:5708590392776 Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com