




PRESSALIT®

Product Item number	Shower mattress for MSCT 1 and MSCT 4 R844721302
Colour	<p>VVS: 777542822 EAN: 5708590373928 HMI: 128434</p> <p>Pewter grey</p> 
Product description	<p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p> <p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p>
Materials	<p>The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.</p>
Warranty	<p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.</p>



PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844721302
Cleaning	<p>MSCT mattresses</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS:</p> <p>Sprit (husholdningssprit)</p> <p>Tvättssprit 70%</p> <p>Safe Zone Plus</p> <p>Selgiene Extreme</p> <p>Microsan</p> <p>Accel PREvention Wipes Lingettes</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®