



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Toilet brush for wall mounting, w/ bowl, brushed steel PY30810
<b>Product description</b>	VVS: 778499115 NRF: 1438655 RSK: 8920393 EAN: 5708590392356  Toilet brush in solid stainless steel for wall mounting. Concealed anti-theft protection with 2-point rosette. Incl. mounting kit.
<b>Materials</b>	European stainless steel (AISI 303)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)