




PRESSALIT®

Product	Wall-hung toilet 700 mm
Item number	R2099000
Colour	VVS:612928000 NRF:6191026 RSK:7822154 EAN:5708590345284 HMI: 105700 White
Product description	Toilet, wall-hung, distance from wall to the front: 700 mm. Excl. mounting kit. Toilet, wall-hung, distance from wall to the front: 700 mm. Excl. mounting kit.
Product capacity	Max load: 400 kg.
Materials	Sanitary porcelain
Ergonomy	700 mm is recommended for wheelchair users as it gives more space for a carer. Wall-hung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



pressalit.com