




PRESSALIT®

Product Item number	Pressalit Choice Toilet brush for wall with glass bowl, brushed steel Q30800
Product description	<p>VVS:778497100 NRF:6190958 RSK:8714261 EAN:5708590337043</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant toilet brush for wall mounting. Opens up extra space on the floor.</p>
Materials	Stainless steel AISI 304, glass, mounting knot in brass
Warranty	A 2-year warranty is provided for material and manufacturing defects.
Cleaning	<p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®