




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, polished steel PY10725
Product description	<p>VVS: 774565640 NRF: 1438616 RSK: 8920383 EAN: 5708590392271</p>  <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com