



PRESSALIT®

Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008

Wall hung toilet with circular flush, incl. toilet seat
Toilet seat with soft close and antibacterial effect, incl. Solid Pro hinge

Item number




CP1008011CF
VVS: 613326300
NRF: 1462950
RSK: 7822114
EAN: 5708590396880

The Pressalit Trend Combi Pack with Projecta D Solid Pro seat is designed for maximum durability and stability. Featuring the innovative rimless Trend bowl with Circular Flush for efficient, splash-free performance, it's the ideal solution for high-demand environments. The Solid Pro seat's robust D-shaped design and reinforced hinges withstand heavy, frequent use and tested with up to 60,000 sideways load cycles. An antibacterial additive helps maintain a hygienic surface and offers uncompromising quality. This combination is perfect for public buildings, institutions, or other settings where strength, stability, and hygiene are paramount, offering a reliable and long-lasting toilet solution.



Colour

White





Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008

Wall hung toilet with circular flush, incl. toilet seat
seal with soft close and antibacterial effect, incl. Solid Pro hinge

Materials**TOILET SEAT:**

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.

ANTIBACTERIAL:

This toilet seat features an antibacterial additive with biocidal silver chloride, which helps prevent bacterial growth.

BUFFERS:

EVA (copolymer made of ethylene and vinyl acetate).

DAMPER KIT (only relevant for toilet seats with soft close):

Hydraulic damper with stainless steel damper casing and thermoplastic parts.

MOUNT/HINGES:

Stainless steel model W. no. 1.4301/M, thermoplastic semi-crystalline polyamides with partially aromatic content. Aluminium profile EN AW-6063-TS.

Product capacity

Load - ring seat 240 kg

Weight

31,508 kg

Carton size

540x410x510 mm

Gross weight

34,223 kg

Weight per crate

34,223 kg

Qty per crate

1

Crates per pallet

5

Qty per pallet

5

Warranty

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:
- supply spare parts, or
- replace the product with an identical or a similar model.
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.



PRESSALIT®

Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008
Wall hung toilet with circular flush, incl. toilet seat
seat with soft close and antibacterial effect, incl. Solid Pro hinge

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®