




# PRESSALIT®

|                            |   |
|----------------------------|---|
| <b>Product</b>             | Toilet paper holder for TMA fold down grab bars   |
| <b>Item number</b>         | R1175000  |
| <b>Color</b>               | <p>VVS:615723950<br/>NRF:6191655<br/>EAN:5708590211725<br/>HMI: 35271<br/>NBS:N13:436<br/>DB: 2076625</p> <p>White</p>  |
| <b>Product description</b> | <p>Paper holder for TMA arm supports.<br/>Fixed, can be mounted on the left or right handle.</p> <p>Paper holder for TMA arm supports.<br/>Fixed, can be mounted on the left or right handle.</p> |
| <b>Adjustability</b>       | Fixed mounting  |
| <b>Materials</b>           | Powder-coated stainless steel.  |
| <b>Ergonomy</b>            | Can be mounted at the right or the left handle.   |
| <b>Warranty</b>            | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.   |



[pressalit.com](http://pressalit.com)

## PRESSALIT®