




PRESSALIT®

Product	Toilet paper holder for TMA fold down grab bars
Item number	R1175000
Color	VVS:615723950 NRF:6191655 RSK:8920763 EAN:5708590211725 HMI: 35271 NBS:N13:436 DB: 2076625 White
Product description	Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle. Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel.
Ergonomy	Can be mounted at the right or the left handle.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



pressalit.com

PRESSALIT®