



United Kingdom Conformity Assessment

United Kingdom RP:
Pressalit Ltd.
100, Longwater Avenue
Green Park
Reading RG2 6GP
United Kingdom
T: 01180923 2696
uk@pressalit.com

This product is produced in compliance with all relevant requirements of:
- Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
- Council Directive 2011/65/EC on RoHS

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019
ISO 13485:2016

The products are manufactured in accordance with the following standards:

ISO 17966:2016

Ry, 23. juni 2026

Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S
Responsible

Manufacturer: Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

VALUE	
REF R102685000	EAN (GTIN): 5708590391069
	Nomenclature (GMDN): 35585
	Nomenclature (EMDN): Y181809
Name	
EN	VALUE support arm with mounting profile, 850 mm, incl. operating panel for electrical toilet lifter, white with grey end cap. Height adjustable, left hand operated, for fitting to TL5 (right)
Intended purpose	
EN	The support arm is used as a suport for a user, when he or she sits down or gets up from a toilet. Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.