



United Kingdom Conformity Assessment

United Kingdom RP:
 Pressalit Ltd.
 100, Longwater Avenue
 Green Park
 Reading RG2 6GP
 United Kingdom
 T: 01180923 2696
 uk@pressalit.com

This product is produced in compliance with all relevant requirements of:
 - Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
 For the UK:
 - Medical Device Regulation (MDR) 2002 on medical devices

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019
 ISO 13485:2016

The products are manufactured in accordance with the following standards:

ISO 17966:2016

Ry, 28. juli 2026

Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S
 Responsible

Manufacturer: Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

PLUS

REF

R4580112

EAN (GTIN): 5708590309415

Nomenclature (GMDN): 63052

Nomenclature (EMDN): Y093399

Name

EN PLUS wash basin bracket, track-mounted, with anthracite cover plate. Manually height adjustable with tilt lock, and manually sideways adjustable

Intended purpose

EN A wash basin bracket is intended for optimum ergonomics at a wash basin, for both seated and standing users with disabilities.
 Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.