



PRESSALIT®

Product	Mounting kit V8663 (8 pcs), for aerated concrete
Item number	V8663
Product description	EAN:5708590351315
	Mounting kit, 8 lag bolts Ø10 x 70, plugs GB14, with grey covers, for aerated concrete.
	Used for: SCT 3000 shower change tables: R8435313, R8435318 Mounting kit, 8 lag bolts Ø10 x 70, plugs GB14, with grey covers, for aerated concrete
Materials	Stainless steel and plastics

pressalit.com