




PRESSALIT®

| | |
|----------------------------|---|
| Product | Shower curtain 1500 mm x 2000 mm, 2 pcs |
| Item number | RT202000 |
| Colour | VVS:777471110 NRF:6311239 EAN:5708590170060 HMI: 42202 NBS:N13:375 White |
| Product description |  Shower curtain, 1500 mm x 2000 mm, 2 singles Shower curtain, 1500 mm x 2000 mm, 2 singles |
| Materials | Polyester cloth with zinc metal eyelet nickelfree and rust free at top, with PE reinforcement lining inside, rubber band at bottom curtain. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. |
| Cleaning | The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®