




PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser with sensor, 850 ml, matt black PY40715
Product description	<p>VVS: 774565231 NRF: 1438640 RSK: 8920398 EAN: 5708590392752</p>  <p>Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.</p>
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com