



PRESSALIT®

Product	Mounting kit V8663 (8 pcs), for aerated concrete
Item number	V8663
Product description	<p>EAN:5708590351315</p> <p>Mounting kit, 8 lag bolts Ø10 x 70, plugs GB14, with grey covers, for aerated concrete.</p> <p>Used for: SCT 3000 shower change tables: R8435313, R8435318</p> <p>Mounting kit, 8 lag bolts Ø10 x 70, plugs GB14, with grey covers, for aerated concrete</p>
Materials	Stainless steel and plastics

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®