



PRESSALIT®

Product Item number	Mounting kit V8652 (2 pcs), for aerated concrete V8652
Product description	EAN:5708590366005 Mounting kit, 2 lag bolts Ø10 x 100, plugs GB14, for aerated concrete. Used for: MATRIX SMALL sinks: R2020/R2022, R2021/R2023 MATRIX CURVE sinks: R2050, R2051, R2052, R2053 Mounting kit, 2 lag bolts Ø10 x 100, plugs GB14, for aerated concrete
Materials	Stainless steel and plastics

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®