




PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717304
	VVS: 777542811 EAN: 5708590373911 HMI: 128437
Color	Pearl gray
Product description	Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 26" x 65.4". Incl. 63" hose. Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 26" x 65.4". Incl. 63" hose.
	
Materials	The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.
Cleaning	This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Sprit (husholdningssprit) Tvättsprit 70% Safe Zone Plus Selgiene Extreme Microsan Accel PREvention Wipes Lingettes The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	Shower mattress for MSCT 1 and MSCT 4
Item number	R844717304

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

8. juli 2026