




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 large roll of toilet paper, brushed steel PX30225
Product description	<p>VVS: 776130115 NRF: 1438649 RSK: 8920367 EAN: 5708590392677</p>  <p>Toilet paper holder in large steel housing with space for rolls of max. Ø 280 x 100 mm. For wall mounting. The holder is lockable and includes mounting kit.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

PRESSALIT®