



PRESSALIT®


Raised buffer with screw 100 mm

Item number

T9804
EAN:5708590183930

Description

Raised buffer with screw 100 mm
Raised buffer with screw 100 mm



Product capacity

Weight

0,1 kg

Gross weight

0,1 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

pressalit.com

PRESSALIT®