



PRESSALIT®

Product Item number	Pressalit Choice Towel hook, pin, 2 pcs., polished steel Q10150
Product description	VVS:775295240 NRF:6190928 RSK:8754882 EAN:5708590336862 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom. Use the towel hook in the bathroom and in the kitchen. That way, you can let the style be consistent throughout the house.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®