




PRESSALIT®

| | |
|----------------------------|---|
| Product | Paper holder for TMA arm supports |
| Item number | R1175000 |
| Colour | VVS:615723950 NRF:6191655 EAN:5708590211725 HMI: 35271 NBS:N13:436 White |
| Product description |  <p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p> <p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p> |
| Adjustability | Fixed mounting |
| Materials | Powder-coated stainless steel. |
| Ergonomy | Can be mounted at the left or right handle. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. |

pressalit.com

PRESSALIT®