




<p>Product Item number</p>	<p>TMA 2 arm support for mounting on toilet, wider fit R1191112000</p>
<p>Color Product description</p>	<p>VVS:615724100 NRF:6190325 EAN:5708590376240 HMI: 130727 DB: 2510918 FDA:ILS/D135866 White</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 2 arm supports for mounting on toilet, incl. BL1999 hinge. Dark gray handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R37 or R38.</p>
<p>Adjustability</p>	<p>Fixed mounting</p>
<p>Materials Ergonomy Features</p>	<p>Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i></p> <p>The arm supports provide firm support and the molded handles are comfortable to hold on to.</p> <p>The arm supports can be folded individually.</p>
<p>Warranty Cleaning</p>	<p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p> <p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>



PRESSALIT®

Product	TMA 2 arm support for mounting on toilet, wider fit
Item number	R1191112000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®