




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Towel holder, single, 60 x 8,5 cm, brushed steel Q30500
<b>Product description</b>	VVS:775711100 NRF:6190955 RSK:8754889 EAN:5708590337029  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p>
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)