



PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser, 950 ml, brushed steel PY30730
Product description	VVS:774565515 NRF:1438617 RSK:8920391 EAN:5708590392332 Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com