



<b>Product Item number</b>	Pressalit Public Classic Paper dispenser, manual, 500 sheets, polished steel PY10110
<b>Product description</b>	VVS:776135540 NRF:1438633 RSK:8920380 EAN:5708590392240  Stainless steel paper dispenser for wall mounting. For approx. 500 paper napkins. Key-operated opening. Incl. mounting kit.  Stainless steel housing. Fills at the top and has an inspection opening to indicate fill level.  Recommended paper size (folded): 240-250 mm x 100-115 mm.
<b>Materials</b>	European stainless steel (AISI 304)
<b>Warranty Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.